



F A R M - T E C

Kett TQ-100 Rice Crusher



If you dry large chunks of samples with an Infrared Moisture Determination Balance, not only will it take a considerable amount of time to heat to the centre of the samples, but the surface of the sample might scorch, making it impossible to get an accurate measurement. What would help you in this case is the TQ-100. The roll-mill was made to crush samples for use in moisture content measuring, and the unit is small enough to use at your own desk. In addition to rice or other grains, this unit is ideal for crushing various food stuffs, aggregates, coal, etc.

Specifications	TQ-100
Sample Weight	Roughly 5g
Grain Size	20-30 mesh (rice)
Dimensions & Weight	80(W)X70(D)X180(H)mm, 1.2kg

Farm-Tec
Unit 1
Carr End Farm
Glaisdale
WHITBY
North Yorkshire
YO21 2QH

Tel: 01947 897995
Fax: 01947 821372

Email:
infofarmtec@aol.com

Website:
www.farm-tec.co.uk